


Product Datasheet

Mouse Cathepsin D ELISA Kit (orb138062)

Description	Mouse Cathepsin D ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I21-L410
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Clonality	Polyclonal
MW	44954 MW
Uniprot ID	Q3UCD9
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 200-300 ng/ml.
Expiration Date	Please enquire.



Mouse
Cathepsin D
PicoKine
ELISA Kit
sta...