


Product Datasheet

Human CA9/Carbonic anhydrase 9 ELISA Kit (orb138063)

Description	Human CA9/Carbonic anhydrase 9 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: P59-D414
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	49698 MW
Uniprot ID	Q16790
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-259 pg/ml.
Expiration Date	Please enquire.



Human CA9
PicoKine
ELISA Kit
standard
cu...