


Product Datasheet

Mouse TNFRSF7/CD27 ELISA Kit (orb138064)

Description	Mouse TNFRSF7/CD27 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T21-R182
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	28164 MW
Uniprot ID	P41272
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
TNFRSF7/CD27
PicoKine ELISA
Kit st...