


Product Datasheet

Mouse TNFSF4/OX40L ELISA Kit (orb138071)

Description	Mouse TNFSF4/OX40L ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S51-L198
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22255 MW
Uniprot ID	P43488
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Mouse
TNFSF4/OX40L
PicoKine ELISA
Kit st...