


Product Datasheet

Human EPCR ELISA Kit (orb138076)

Description	Human EPCR ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S18-S210
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	26671 MW
Uniprot ID	Q9UNN8
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is around 160 ng/ml.
Expiration Date	Please enquire.



Human EPCR
PicoKine
ELISA Kit
standard c...