


Product Datasheet

Human ANGPTL4 / Angiopoietin Like 4 / ARP4 ELISA Kit (orb138079)

Description	Human ANGPTL4 / Angiopoietin Like 4 / ARP4 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L165-S406
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	45214 MW
Uniprot ID	Q9BY76
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1.2 ng/ml.
Expiration Date	Please enquire.



Human
ANGPTL4
PicoKine
ELISA Kit
standar...