


Product Datasheet

Human IGFBP6 ELISA Kit (orb138081)

Description	Human IGFBP6 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A25-G240
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	25322 MW
Uniprot ID	P24592
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is 197-438 ng/ml.
Expiration Date	Please enquire.



Human IGFBP6 ELISA Kit standard curve