


Product Datasheet

Human IL-29 ELISA Kit (orb138083)

Description	Human IL-29 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G20-T200
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	21898 MW
Uniprot ID	Q8IU54
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-29
PicoKine
ELISA Kit
standard ...