


Product Datasheet

Human complement C5a ELISA Kit (orb138085)

Description	Human complement C5a ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T678-R751
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	188305 MW
Uniprot ID	P01031
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 80 ng/ml.
Expiration Date	Please enquire.



Human
complement
C5a
PicoKine
ELISA Kit ...