


Product Datasheet

Human LYVE-1 ELISA Kit (orb138086)

Description	Human LYVE-1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S24-T238
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and cell lysates.
MW	35213 MW
Uniprot ID	Q9Y5Y7
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum is around 50ng/ml
Expiration Date	Please enquire.



Human
LYVE-1
PicoKine
ELISA Kit
standard...