


Product Datasheet

Human TNFRSF13B/TACI ELISA Kit (orb138089)

Description	Human TNFRSF13B/TACI ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S2-T166
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	31816 MW
Uniprot ID	O14836
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 300 pg/ml.
Expiration Date	Please enquire.



Human
TNFRSF13B/TACI
PicoKine ELISA
Kit ...