


Product Datasheet

Mouse Complement C5a ELISA Kit (orb138090)

Description	Mouse Complement C5a ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N679-R755
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	188878 MW
Uniprot ID	P06684
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 10 ng/ml.
Expiration Date	Please enquire.



Mouse
Complement
C5a
PicoKine
ELISA Kit ...