


Product Datasheet

Human TNFRSF9/4-1BB ELISA Kit (orb138091)

Description	Human TNFRSF9/4-1BB ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L24-Q186
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27899 MW
Uniprot ID	Q07011
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TNFRSF9/4-
1BB
PicoKine
ELISA Kit s...