


Product Datasheet

Human TNFSF14/LIGHT ELISA Kit (orb138092)

Description	Human TNFSF14/LIGHT ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D74-V240
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	26350 MW
Uniprot ID	O43557
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-259 pg/ml.
Expiration Date	Please enquire.



Human
TNFSF14/LIGHT
PicoKine ELISA
Kit s...