


Product Datasheet

Human IGFBP7/Igfbp Rp1 ELISA Kit (orb138093)

Description	Human IGFBP7/Igfbp Rp1 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D30-L282
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	29130 MW
Uniprot ID	Q16270
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 50 ng/ml.
Expiration Date	Please enquire.



Human
IGFBP7
PicoKine
ELISA Kit
standard...