


Product Datasheet

Human CD200 ELISA Kit (orb138094)

Description	Human CD200 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q31-G232
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	31264 MW
Uniprot ID	P41217
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CD200
PicoKine
ELISA Kit
standard ...