


Product Datasheet

Human ICAM-2 ELISA Kit (orb138100)

Description	Human ICAM-2 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K25-Q223
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	30654 MW
Uniprot ID	P13598
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 10.4-41.2 ng/ml.
Expiration Date	Please enquire.



Human
ICAM-2
PicoKine
ELISA Kit
standard...