


Product Datasheet

Rat FASL ELISA Kit (orb138105)

Description	Rat FASL ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q127-L278
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	31140 MW
Uniprot ID	P36940
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat FASL
PicoKine
ELISA Kit
standard
cur...