


Product Datasheet

Human Galectin-9 / Gal 9 ELISA Kit (orb138106)

Description	Human Galectin-9 / Gal 9 ELISA Kit
Reactivity	Human
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-T355
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39518 MW
Uniprot ID	000182
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 2.48-12.48 ng/ml.
Expiration Date	Please enquire.



Human
Galectin-9
PicoKine
ELISA Kit
stan...