



Product Datasheet

Aurora B antibody (FITC) (orb13817)


Description	Aurora B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Aurora B
Target	Aurora B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96GD4
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(human) lys...



Western blot analysis of Bone (Mouse) Ly...



Flow cytometric analysis of A431 cells U...