

Product Datasheet

GUCY1A2 peptide (orb13867)

Description	GUCY1A2 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	GUCY1A2
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb389331 (GUCY1A2 antibody)
Expiration Date	6 months from date of receipt.