

Product Datasheet

ER alpha (phospho-Ser167) antibody (FITC) (orb13967)

Description	ER alpha (phospho-Ser167) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Estrogen Receptor alpha around the phosphorylation site of Ser167
Target	Phospho-Estrogen Receptor alpha (Ser167)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	Flow-Cyt=0.2µg /test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.