

Product Datasheet

MMP9 antibody (orb140003)

Description	Mouse Monoclonal to MMP9.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MMP9
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: MMP9. Gene-IDs: 4318. Accessions: NP_004985.2, P14780.
Isotype	IgM
Clonality	Monoclonal
MW	78458
Uniprot ID	P14780
NCBI	NP_004985.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

