

Product Datasheet

TSC1 antibody (orb140005)

Description	Mouse Monoclonal to TSC1.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	TSC1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: TSC1. Gene-IDs: 7248. Accessions: NP_000359.1, NP_001155898.1, Q92574.
Isotype	IgG2b
Clonality	Monoclonal
MW	129767
Uniprot ID	Q92574
NCBI	NP_001155898.1 , NP_001155899.1 , NP_000359.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

