

Product Datasheet

TLE1 antibody (orb140007)

Description	Mouse Monoclonal to TLE1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TLE1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: TLE1. Gene-IDs: 7088. Accessions: Q04724, NP_005068.2.
Isotype	IgG1
Clonality	Monoclonal
MW	83201
Uniprot ID	Q04724
NCBI	NP_005068.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

