

Product Datasheet

INS antibody (orb140009)

Description	Mouse Monoclonal to INS.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	INS
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: INS. Gene-IDs: 3630. Accessions: NP_000198.1, NP_001172027.1, NP_001172026.1.
Isotype	IgM
Clonality	Monoclonal
MW	11981
Uniprot ID	P01308
NCBI	NP_001172026.1, NP_001172027.1, NP_000198.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

