

Product Datasheet

RPS6KA3 antibody (orb140019)

Description	Mouse Monoclonal to RPS6KA3.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RPS6KA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: RPS6KA3. Gene-IDs: 6197. Accessions: P51812, NP_004577.1.
Isotype	IgG1
Clonality	Monoclonal
MW	83736
Uniprot ID	P51812
NCBI	NP_004577.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

