

Product Datasheet

MUC1 antibody (orb140020)

Description	Mouse Monoclonal to MUC1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MUC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: MUC1. Gene-IDs: 4582. Accessions: NP_001018017.1, P15941, NP_001018016.1.
Isotype	IgG2b
Clonality	Monoclonal
MW	122102
Uniprot ID	P15941
NCBI	NP_001191218.1, NP_001191217.1, NP_001037858.1, NP_001018016.1, NP_001018017.1, NP_00, NP_001191216.1, NP_001037856.1, NP_001191215.1, NP_001037855.1, NP_001191219.1, NP_001037857.1, NP_001191220.1, NP_001191214.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

