

Product Datasheet

TLE4 antibody (orb140027)

Description	Mouse Monoclonal to TLE4.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TLE4
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: TLE4. Gene-IDs: 7091. Accessions: Q04727, NP_008936.2.
Isotype	IgM
Clonality	Monoclonal
MW	83755
Uniprot ID	Q04727
NCBI	NP_008936.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

