

Product Datasheet

HNMT antibody (orb140032)

Description	Mouse Monoclonal to HNMT.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	HNMT
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: HNMT. Gene-IDs: 3176. Accessions: NP_008826.1, P50135.
Isotype	IgG1
Clonality	Monoclonal
MW	33295
Uniprot ID	P50135
NCBI	NP_001019245.1, NP_001019246.1, NP_008826.1
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

