

Product Datasheet

RAB23 antibody (orb140037)

Description	Mouse Monoclonal to RAB23.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant Protein
Target	RAB23
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: RAB23. Gene-IDs: 51715. Accessions: Q9ULC3, NP_057361.3, NP_899050.1.
Isotype	IgG1
Clonality	Monoclonal
MW	26659
Uniprot ID	Q9ULC3
NCBI	NP_001265597.1 , NP_899050.1 , NP_001265595.1 , NP_057361.3 , NP_001265596.1
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

