



Product Datasheet

GP2 antibody (orb140051)

Description	Rabbit Polyclonal to GP2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: GP2. Gene-IDs: 2813. Accessions: NP_001007241.2, NP_001493.2, NP_001007243.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59480
Uniprot ID	P55259
NCBI	NP_001007242.2 , NP_001493.2 , NP_001007243.2 , NP_001007241.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...