



Product Datasheet

LIMD1 antibody (orb140055)

Description	Rabbit Polyclonal to LIMD1.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	LIMD1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: LIMD1. Gene-IDs: 8994. Accessions: Q9UGP4, NP_055055.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72190
Uniprot ID	Q9UGP4
NCBI	NP_055055.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...