



Product Datasheet

DNAJC14 antibody (orb140064)

Description	Rabbit Polyclonal to DNAJC14.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	DNAJC14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: DNAJC14. Gene-IDs: 85406. Accessions: Q6Y2X3, NP_115740.5.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78569
Uniprot ID	Q6Y2X3
NCBI	NP_115740.5
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...