



Product Datasheet

GGT7 antibody (orb140069)

Description	Rabbit Polyclonal to GGT7.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GGT7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: GGT7. Gene-IDs: 2686. Accessions: NP_821158.2, Q9UJ14.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70467
Uniprot ID	Q9UJ14
NCBI	NP_821158.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...