




Product Datasheet

PSORS1C1 antibody (orb140087)

Description	Rabbit Polyclonal to PSORS1C1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PSORS1C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PSORS1C1. Gene-IDs: 170679. Accessions: NP_054787.2, Q9UIG5.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16580
Uniprot ID	Q9UIG5
NCBI	NP_054787.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of K562 cells (...)