



Product Datasheet

CD11c antibody (FITC) (orb14009)

Description	CD11c antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Integrin Alpha X
Target	CD11c
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20702
Dilution Range	Flow-Cyt=1µg/Test,IF=1:200-800
Expiration Date	12 months from date of receipt.



Blank control:
Mouse
spleen.
Primary Ant...



Blank
control:Mouse
spleen.
Primary
Anti...