




Product Datasheet

FGF11 antibody (orb140092)

Description	Rabbit Polyclonal to FGF11.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	FGF11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: FGF11. Gene-IDs: 2256. Accessions: NP_004103.1, Q92914.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25005
Uniprot ID	Q92914
NCBI	NP_004103.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of A2058 cells ...