



Product Datasheet

FOXA2 antibody (orb140093)

Description	Rabbit Polyclonal to FOXA2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	FOXA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: FOXA2. Gene-IDs: 3170. Accessions: Q9Y261, NP_068556.2, NP_710141.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48306
Uniprot ID	Q9Y261
NCBI	NP_710141.1 , NP_068556.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...