



Product Datasheet

NDST1 antibody (orb140096)

Description	Rabbit Polyclonal to NDST1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NDST1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: NDST1. Gene-IDs: 3340. Accessions: P52848, NP_001534.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	100868
Uniprot ID	P52848
NCBI	NP_001534.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of NDST1 antibody ...