



Product Datasheet

GPHB5 antibody (orb140101)

Description	Rabbit Polyclonal to GPHB5.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GPHB5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: GPHB5. Gene-IDs: 122876. Accessions: Q86YW7, NP_660154.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14232
Uniprot ID	Q86YW7
NCBI	NP_660154.3
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...