



Product Datasheet

ERO1L antibody (orb140121)

Description	Rabbit Polyclonal to ERO1L.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ERO1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ERO1L. Gene-IDs: 30001. Accessions: NP_055399.1, Q96HE7.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54393
Uniprot ID	Q96HE7
NCBI	NP_055399.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...