



Product Datasheet

WIPI2 antibody (orb140123)

Description	Rabbit Polyclonal to WIPI2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	WIPI2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: WIPI2. Gene-IDs: 26100. Accessions: NP_001028690.1, NP_057087.2, NP_001028691.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49408
Uniprot ID	Q9Y4P8
NCBI	NP_001265228.1 , NP_001028692.1 , NP_001028690.1 , NP_057087.2 , NP_001028691.1 , NP_056425.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...