



Product Datasheet

SGEF antibody (orb140129)

Description	Rabbit Polyclonal to SGEF.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ARHGEF26
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ARHGEF26. Gene-IDs: 26084. Accessions: NP_056410.3, Q96DR7.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97346
Uniprot ID	Q96DR7
NCBI	NP_001238892.1, NP_001238891.1, NP_056410.3
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...