



Product Datasheet

ZNF227 antibody (orb140130)

Description	Rabbit Polyclonal to ZNF227.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ZNF227
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ZNF227. Gene-IDs: 7770. Accessions: Q86WZ6, NP_872296.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	92033
Uniprot ID	Q86WZ6
NCBI	NP_872296.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...