




Product Datasheet

PI15 antibody (orb140132)

Description	Rabbit Polyclonal to PI15.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PI15
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PI15. Gene-IDs: 51050. Accessions: O43692, NP_056970.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	29065
Uniprot ID	O43692
NCBI	NP_056970.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of K562 cells (...)