




Product Datasheet

ELAVL3 antibody (orb140134)

Description	Rabbit Polyclonal to ELAVL3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ELAVL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ELAVL3. Gene-IDs: 1995. Accessions: Q14576, NP_115657.2, NP_001411.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39547
Uniprot ID	Q14576
NCBI	NP_115657.2 , NP_001411.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of ELAVL3 antibody...



WB analysis of human brain tissue lysate...