



Product Datasheet

ROBLD3 antibody (orb140135)

Description	Rabbit Polyclonal to ROBLD3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	LAMTOR2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ROBLD3. Gene-IDs: 28956. Accessions: NP_054736.1, Q9Y2Q5.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13508
Uniprot ID	Q9Y2Q5
NCBI	NP_054736.1 , NP_001138736.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...