



Product Datasheet

DHRS7B antibody (orb140136)

Description	Rabbit Polyclonal to DHRS7B.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	DHRS7B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: DHRS7B. Gene-IDs: 25979. Accessions: Q6IAN0, NP_056325.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35119
Uniprot ID	Q6IAN0
NCBI	NP_056325.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...