



Product Datasheet

SLIT1 antibody (orb140141)

Description	Rabbit polyclonal antibody to SLIT3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SLIT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SLIT1. Gene-IDs: 6585. Accessions: O75093, NP_003052.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	167926
Uniprot ID	O75093
NCBI	NP_003052.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...